

## FACS Aria Fusion Filter Setup

Channel		Bandpass	Fluorochromes
<b>VIOLET LASER (405nm)</b>			
C		450/50	Pacific Blue, CFP, e450, BFP, BV421, DAPI, Zombie Violet, Alexa Fluor 405, mCerulean, mTurquoise, Cascade Blue, DyeCycle Violet, Hoechst L/D
B		525/50	BV510, L/D AQUA, Saffire, Pacific Orange, AmCyan
A		610/20	BV605, BV570, Zombie Yellow, Superbright 600
<b>YELLOW LASER (561 nm)</b>			
D		582/15	PE, DsRed, mOrange, mTomato, CY3, PI, Alexa Fluor 555, RFP
C		610/20	mCherry, PE-CF594, PE-Texas Red, mRuby, mScarlet
B		670/14	PE-Cy5, PE-CY5.5, mKate, mPlum
A		780/60	PE-Cy7
<b>RED LASER (638 nm)</b>			
B		670/30	APC, Alexa Fluor 647, Sytox Red, iRFP
A		780/60	APC-Cy7, APC-e780, APC-H7, AF 750
<b>BLUE LASER (488 nm)</b>			
B		530/30	GFP, FITC, Alexa Fluor 488, YFP, Neon Green, Venus, mCitrine, CFSC
A		695/40	7AAD, PerCP-Cy5.5, PE Cy5.5