

Miltenyi VYB Filter Setup

Channel		Bandpass	Fluorochromes
VIOLET LASER (405nm)			
V1		450/40	Pacific Blue, CFP, e450, BFP, BV421, DAPI, Zombie Violet, AF405, mCerulean, mTurquoise, Cascade Blue, DyeCycle Violet, Hoechst L/D, AmCyan
V2		525/50	BV510, CFP/YFP Fret, L/D AQUA, Saffire, Pacific Orange, VioGreen, Sytox Blue
YELLOW LASER (561 nm)			
Y1		586/15	PE, mTomato
Y2		610/20	mCherry, Texas Red, DsRed
Y3		661/20	APC, PECy5, mKate, PerCPCY5.5
Y4		750 LP	APC-Cy7, APC-e780, PE-CY7, APC-H7
BLUE LASER (488 nm)			
B1		525/50	GFP, FITC, YFP
B2		614/50	LSS-mKate, YFP, PI