

FACS Aria Fusion Filter Setup

Channel	Name	Bandpass	Fluorochromes
VIOLET LASER (405nm)			
C	V450	450/50	Pacific Blue, CFP, e450, BFP, BV421, DAPI, Zombie Violet, Alexa Fluor 405, mCerulean, mTurquoise, Cascade Blue, DyeCycle Violet, Hoechst L/D
B	V525	525/50	BV510, L/D AQUA, Saffire, Pacific Orange, AmCyan
A	V510	610/20	BV605, BV570, Zombie Yellow, Superbright 600
YELLOW LASER (561 nm)			
D	Y582	582/15	PE, DsRed, mOrange, mTomato, CY3, PI, Alexa Fluor 555, RFP
C	Y610	610/20	mCherry, PE-CF594, PE-Texas Red, mRuby, mScarlet
B	Y670	670/40	PE-Cy5, PE-CY5.5, mKate, mPlum
A	Y780	780/60	PE-Cy7
RED LASER (638 nm)			
B	R670	670/30	APC, Alexa Fluor 647, Sytox Red, iRFP
A	R780	780/60	APC-Cy7, APC-e780, APC-H7, AF 750
BLUE LASER (488 nm)			
B	B530	530/30	GFP, FITC, Alexa Fluor 488, YFP, Neon Green, Venus, mCitrine, CFSC
A	B695	695/40	7AAD, PerCP-Cy5.5, PE Cy5.5